





# REPLACEMENT BP CUFF REFERENCE GUIDE

	CSM and ProBP 3400			CVSM and SPOT LXI			Atlas, VSM and Spot Devices (not LXI)		
Connector	No Connector (Raw)			2-tube Locking Type			1-tube Screw-Type		
									
Type	Durable	Disposable	Disposable	Durable	Disposable	Disposable	Durable	Disposable	Disposable
Material	Nylon	Soft	Vinyl	Nylon	Soft	Vinyl	Nylon	Soft	Vinyl
Unit of Measure	Each	20 Pack	20 Pack	Each	20 Pack	20 Pack	Each	20 Pack	20 Pack
Sm. Infant	Reuse-06	Soft-06	Vinyl-06	N/A	N/A	N/A	Reuse-06-1SC	Soft-06-1SC	Vinyl-06-1SC
Infant	Reuse-07	Soft-07	Vinyl-07	Reuse-07-2MQ	Soft-07-2MQ	Vinyl-07-2MQ	Reuse-07-1SC	Soft-07-1SC	Vinyl-07-1SC
Sm. Child	Reuse-08	Soft-08	Vinyl-08	Reuse-08-2MQ	Soft-08-2MQ	Vinyl-08-2MQ	Reuse-08-1SC	Soft-08-1SC	Vinyl-08-1SC
Child	Reuse-09	Soft-09	Vinyl-09	Reuse-09-2MQ	Soft-09-2MQ	Vinyl-09-2MQ	Reuse-09-1SC	Soft-09-1SC	Vinyl-09-1SC
Sm. Adult	Reuse-10	Soft-10	Vinyl-10	Reuse-10-2MQ	Soft-10-2MQ	Vinyl-10-2MQ	Reuse-10-1SC	Soft-10-1SC	Vinyl-10-1SC
Adult	Reuse-11	Soft-11	Vinyl-11	Reuse-11-2MQ	Soft-11-2MQ	Vinyl-11-2MQ	Reuse-11-1SC	Soft-11-1SC	Vinyl-11-1SC
Adult Long	Reuse-11L	Soft-11L	Vinyl-11L	Reuse-11L-2MQ	Soft-11L-2MQ	Vinyl-11L-2MQ	Reuse-11L-1SC	Soft-11L-1SC	Vinyl-11L-1SC
Lg. Adult	Reuse-12	Soft-12	Vinyl-12	Reuse-12-2MQ	Soft-12-2MQ	Vinyl-12-2MQ	Reuse-12-1SC	Soft-12-1SC	Vinyl-12-1SC
Lg. Adult Long	Reuse-12L	Soft-12L	Vinyl-12L	Reuse-12L-2MQ	Soft-12L-2MQ	Vinyl-12L-2MQ	Reuse-12L-1SC	Soft-12L-1SC	Vinyl-12L-1SC
Thigh	Reuse-13	Soft-13	Vinyl-13	Reuse-13-2MQ	Soft-13-2MQ	Vinyl-13-2MQ	Reuse-13-1SC	Soft-13-1SC	Vinyl-13-1SC

	Propaq and PIC Monitors			Philips, Hewlett Packard (HP), Spacelabs, Datascope and Mindray Monitors			Critikon Dinamap Monitors		
Connector	1-tube Locking Type			1-tube Bayonet (aka HP, Rectus)			2-tube Screw-Type		
									
Type	Durable	Disposable	Disposable	Durable	Disposable	Disposable	Durable	Disposable	Disposable
Material	Nylon	Soft	Vinyl	Nylon	Soft	Vinyl	Nylon	Soft	Vinyl
Unit of Measure	Each	20 Pack	20 Pack	Each	20 Pack	20 Pack	Each	20 Pack	20 Pack
Sm. Infant	Reuse-06-1MQ	Soft-06-1MQ	Vinyl-06-1MQ	Reuse-06-1HP	Soft-06-1HP	Vinyl-06-1HP	N/A	N/A	N/A
Infant	Reuse-07-1MQ	Soft-07-1MQ	Vinyl-07-1MQ	Reuse-07-1HP	Soft-07-1HP	Vinyl-07-1HP	Reuse-07-2SC	Soft-07-2SC	Vinyl-07-2SC
Sm. Child	Reuse-08-1MQ	Soft-08-1MQ	Vinyl-08-1MQ	Reuse-08-1HP	Soft-08-1HP	Vinyl-08-1HP	Reuse-08-2SC	Soft-08-2SC	Vinyl-08-2SC
Child	Reuse-09-1MQ	Soft-09-1MQ	Vinyl-09-1MQ	Reuse-09-1HP	Soft-09-1HP	Vinyl-09-1HP	Reuse-09-2SC	Soft-09-2SC	Vinyl-09-2SC
Sm. Adult	Reuse-10-1MQ	Soft-10-1MQ	Vinyl-10-1MQ	Reuse-10-1HP	Soft-10-1HP	Vinyl-10-1HP	Reuse-10-2SC	Soft-10-2SC	Vinyl-10-2SC
Adult	Reuse-11-1MQ	Soft-11-1MQ	Vinyl-11-1MQ	Reuse-11-1HP	Soft-11-1HP	Vinyl-11-1HP	Reuse-11-2SC	Soft-11-2SC	Vinyl-11-2SC
Adult Long	Reuse-11L-1MQ	Soft-11L-1MQ	Vinyl-11L-1MQ	Reuse-11L-1HP	Soft-11L-1HP	Vinyl-11L-1HP	Reuse-11L-2SC	Soft-11L-2SC	Vinyl-11L-2SC
Lg. Adult	Reuse-12-1MQ	Soft-12-1MQ	Vinyl-12-1MQ	Reuse-12-1HP	Soft-12-1HP	Vinyl-12-1HP	Reuse-12-2SC	Soft-12-2SC	Vinyl-12-2SC
Lg. Adult Long	Reuse-12L-1MQ	Soft-12L-1MQ	Vinyl-12L-1MQ	Reuse-12L-1HP	Soft-12L-1HP	Vinyl-12L-1HP	Reuse-12L-2SC	Soft-12L-2SC	Vinyl-12L-2SC
Thigh	Reuse-13-1MQ	Soft-13-1MQ	Vinyl-13-1MQ	Reuse-13-1HP	Soft-13-1HP	Vinyl-13-1HP	Reuse-13-2SC	Soft-13-2SC	Vinyl-13-2SC

	GE Marquette Monitors		
Connector	2-tube Locking-Type		
Type	Durable	Disposable	Disposable
Material	Nylon	Soft	Vinyl
Unit of Measure	Each	20 Pack	20 Pack
Sm. Infant	N/A	N/A	N/A
Infant	Reuse-07-2MQ	Soft-07-2MQ	Vinyl-07-2MQ
Sm. Child	Reuse-08-2MQ	Soft-08-2MQ	Vinyl-08-2MQ
Child	Reuse-09-2MQ	Soft-09-2MQ	Vinyl-09-2MQ
Sm. Adult	Reuse-10-2MQ	Soft-10-2MQ	Vinyl-10-2MQ
Adult	Reuse-11-2MQ	Soft-11-2MQ	Vinyl-11-2MQ
Adult Long	Reuse-11L-2MQ	Soft-11L-2MQ	Vinyl-11L-2MQ
Lg. Adult	Reuse-12-2MQ	Soft-12-2MQ	Vinyl-12-2MQ
Lg. Adult Long	Reuse-12L-2MQ	Soft-12L-2MQ	Vinyl-12L-2MQ
Thigh	Reuse-13-2MQ	Soft-13-2MQ	Vinyl-13-2MQ

**WelchAllyn**<sup>®</sup>