












































Lamp Number	Replacement Lamp	Current Instruments	Obsolete Instruments
04700-U		<p>Standard Laryngoscopes</p> <p>680 Series Sizes 0-1 </p> <p>690 Series Sizes 1-2 </p> <p>692 Series Sizes 1-2 </p>	<p>Standard Laryngoscopes</p> <p>634 Series Sizes 0-1 </p> <p>664 Series Size 1 </p> <p>674 Series Size 1 </p> <p>684 Series Sizes 0-1 </p> <p>694 Series Sizes 1-2 </p>
04800-U		<p>Standard Laryngoscopes</p> <p>680 Series Sizes 2-4 </p> <p>690 Series Sizes 3-4 </p> <p>692 Series Sizes 3-4 </p>	<p>Standard Laryngoscopes</p> <p>634 Series Sizes 2-4 </p> <p>640 Series Size 1 </p> <p>640 Series Sizes 2-4 </p> <p>664 Series Sizes 2-4 </p> <p>674 Series Sizes 2-4 </p> <p>684 Series Sizes 2-4 </p> <p>694 Series Sizes 3-4 </p> <p>Transilluminator 40415 </p>
04900-U		<p>Ophthalmoscopes</p> <p> 11720</p> <p> 11730/11735</p>	<p>Ophthalmoscopes</p> <p> 11620</p> <p> 11630</p>
06000-U		<p>Laryngoscope Handles</p> <p> 60813</p> <p> 60814</p> <p> 60713</p> <p> 60815</p>	<p>Laryngoscope Handles</p> <p> 60803</p> <p> 60804</p>
06100-U			<p>Examination Lights</p> <p>48400/48435 </p>
06200-U		<p>AudioScope</p> <p>23300 </p>	<p>AudioScope</p> <p>23000, 23020, 23040 </p>
06300-U			<p>Examination Lights</p> <p>44100 </p> <p>44110 </p>
06400-U			<p>Examination Light</p> <p>44200 </p>